



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Steingrabenstrasse 14  
DE-82211 Herrsching  
Germany  
Delivery Address: Am Anger 5, DE-82237 Woerthsee-Etterschlag

Tel +49 (0) 8153 88702-00  
Fax +49 (0) 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)


## Standard 32Ch actiCAP Electrode Cap for LiveAmp

● Green Holders: Label 1-32, Hardware Channel 1 – 32

● Blue holder: Label & hard-wired Ref

● Black holder: Label & hard-wired Gnd

### Electrode Names and Number Labels



## Details

### Ordering Information

For ordering please give Article Number, Cap Cut, and Size  
(**Example: CLA-32, Caucasian, 58**):

- Article Number: **CLA-32**
- Cap Cut: **Caucasian** or **Asian**
- Size (given in cm head circumference):
  - Adult caps: **54, 56, 58, 60, 62, 64** (average male: 58, average female: 54)
  - Children caps: **50** (5 years), **52** (7 years), **54** (11 years)
  - Infant caps: **36** (newborn), **38, 40** (3 months), **42, 44** (7 month), **46, 48** (3 years)

The catalogue-number comprises the cap as described, serial number, and this document; all packed in a labelled cardboard box. For further information about accessories or consumables please visit our website or contact our local distributor.

### Cap

White Sub-Inion Cap with Pocket for LiveAmp and integrated chin belt,  
Sizes 52 – 64 made from High-Precision Fabric, Sizes 50 and smaller made from High Comfort Fabric  
*Options: Caucasian or Asian, Size*

The holders are mounted according to the figure. The positions include the 19 sites of the “International 10/20-System”. The core positions are extended by further positions of the 10%-System to increase the topographical resolution. Positions TP9/TP10 approximately meet left and right mastoids.

### Accessories

The cap can be fixed either by the integrated chin belt alone, or additionally / alternatively by a chest belt (not included). Normally the decision between these two methods can be made according to what is more comfortable for the subject. Under difficult conditions the chest belt will hold the cap in place a little better.

### Service

In case of any further questions please contact our local distributor or [info@easycap.de](mailto:info@easycap.de). Please have the serial number of your cap at hand (syntax CLAxxxxx, located on the backside label of the cap).